


Medical Gas Valves & the ASHE PDC Summit

For use in medical gas systems, pharmaceuticals and other applications where brazing or soldering is the preferred method of installation, Milwaukee Valve offers a 3-piece bronze full-port ball valve with 6" copper tube extensions brazed to each tailpiece. Because there is no need to disassemble the valve as required when installing other 3-piece valves, contractors can realize considerable labor savings with the [BA350TE](#).



The [BA350TE](#) is available in sizes ranging from 1/2" – 3". For increased installation flexibility, the [BA350TE](#) can be ordered with straight male tube-ends, female expanded ends or with a combination of male and female ends. Manufactured and tested in accordance with MSS SP-110, the valve's full-port design provides high flow with minimal pressure drop, and the 3-piece construction simplifies maintenance and repair.

For more information, look up the [BA350TE](#) at www.MilwaukeeValve.com or contact your sales representative or regional manager. You can also see the [BA350TE](#) at the [ASHE PDC Summit](#) next week in New Orleans. The PDC Summit is the premiere event for health care and hospital facility Planning, Design and Construction decision-makers.



Swing by the New Orleans Ernest Morial Convention Center and visit us at [Booth 911](#). The tradeshow portion of the Summit runs from:

- Monday, March 21
 - 10:30AM – 2:00PM & 4:30PM – 6:00PM
- Tuesday, March 22
 - 10:30AM – 2:00PM



MILWAUKEE VALVE

www.milwaukeevalve.com
262.432.2700



HAMMOND VALVE

www.hammondvalve.com
262.432.2702



www.milwaukeevalve.com
262.432.2700